

How to store and handle your medication

KINIKSA
oneconnect[®]

The Kiniksa OneConnect program provides you with product support for treatment with ARCALYST[®] (rilonacept). Medication storage is an important part of your treatment experience, and as your dedicated Patient Access Leads, we can help you understand how to keep your ARCALYST ready and safe to use.

Storing your ARCALYST

It is important to properly store your ARCALYST at all times.

- ✓ Keep ARCALYST in the carton it comes in to protect it from light.
- ✓ Store ARCALYST in the refrigerator between 36°F and 46°F (2°C and 8°C).
- ✓ ARCALYST may be kept at room temperature after mixing and should be used within 3 hours of mixing. Keep it away from light.



Handling your ARCALYST

In order to prioritize safety, be aware of the expiration date of your ARCALYST, which supplies you are using for each injection, and where those items are to be discarded after use.

- ✓ Refrigerated ARCALYST can be used until the expiration date printed on the vial and carton.
- ✓ Properly dispose of your used vials after injection by placing them in your provided sharps disposal container, including any unused ARCALYST.
- ✓ Keep ARCALYST, injection supplies, and all other medicines out of reach of children.

We're here to help.

If you have any questions, don't hesitate to contact the Kiniksa OneConnect support program team.

1-833-KINIKSA (1-833-546-4572) • Monday-Friday (8 AM-8 PM ET)

Learn more about our personalized treatment support at [KiniksaOneConnect.com](https://www.kiniksaoneconnect.com)



ARCALYST is a registered trademark of Regeneron Pharmaceuticals, Inc.
Kiniksa OneConnect is a registered trademark of Kiniksa Pharmaceuticals.
©2026 Kiniksa Pharmaceuticals. All Rights Reserved. 05/26 ARC-US-00127-23 v3

Arcalyst[®]
(rilonacept) For Injection